

Tracheobronchial Nebulizer

Nebulizer and Kits

SINGLE - USE



MEDICAL EQUIPMENT



Area of application

Nebulizers are used to provide aerosolized humidification or medication for patients to inhale. The product is for single use only.

Product range

Tube Length	Connector (to Oxygen Source)	Mask Size	Description	Reference
Tracheobronchial Nebulizer				
-	-	-	Nebulizer Chamber only (No tubing)	93-770Q
Tracheobronchial Nebulizer Handheld Kits				
2.1m	Universal Connector	-	Hand-held Kit with T-Adapter, Mouthpiece, 15cm reservoir tube, 2.1m (7') Oxygen safety tube, Universal connector	93-771Q
2.1m	Standard Connector	-	Hand-held Kit with T-Adapter, Mouthpiece, 15cm reservoir tube, 2.1m (7') Oxygen safety tube, Standard connector	9-771Q
Tracheobronchial Nebulizer Mask-Kits				
2.1m	Universal Connector	Adult	Nebulizer Kit, Adult aerosol mask (elongated) with elastic strap, 2.1m (7') Oxygen safety tube, Universal connector	93-772Q
2.1m	Standard Connector	Adult	Nebulizer Kit, Adult aerosol mask (elongated) with elastic strap, 2.1m (7') Oxygen safety tube, Universal connector	9-772Q 9-772QV (Vanilla scented mask)
2.1m	Universal Connector	Pediatric	Nebulizer Kit, Pediatric aerosol mask (elongated) with elastic strap, 2.1m (7') Oxygen safety tube, Universal connector	93-773Q
2.1m	Standard Connector	Pediatric	Nebulizer Kit, Pediatric aerosol mask (elongated) with elastic strap, 2.1m (7') Oxygen safety tube, Standard connector	9-773Q 9-773QV (Vanilla scented mask)

Sterilization

Non-sterile, Shelf life: 5 years

Instructions for use

Included? Yes



Tracheobronchial Nebulizer

Nebulizer and Kits

Certifications	
CE-marking	CE 1639
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485 Medical Device Single Audit Program (MDSAP)

Material specifications	
Contains latex?	No
Contains phthalates?	No
Contains Bisphenols?	No
Contains cotton?	No
Contains paper?	Packaging only
Aerosol mask	Polyvinyl chloride
Elastic strap	Polyester (Latex-free)
Aerosol connector	Polypropylene
Nebulizer	K-Resin styrene-butadiene copolymers
Aerosol/reservoir tube	Ethylvinylacetate
Mouthpiece	High Density polyethylene
T-Adaptor	Acrylonitrile-butadiene-styrene
Oxygen tube	Polyvinyl chloride
Connector	Polyvinyl chloride

Packaging information							
REF	Primary package/Polybag	Secondary package/Polybag	Tertiary pack./Shipper			Gross weight kg	Pallet Qty. pcs
	Qty. pcs	Qty. pcs	Qty. pcs	Dimensions mm (LxWxD)	Volume m ³		
93-770Q	1	10	50	300 x 200 x 200	0.021	1.46	8,000
93-771Q	1	10	50	400 x 350 x 320	0.045	4.40	1,800
9-771Q	1	10	50	400 x 350 x 320	0.045	4.40	1,800
93-772Q	1	10	50	400 x 350 x 320	0.045	5.20	1,800
9-772Q	1	10	50	400 x 350 x 320	0.045	5.20	1,800
9-772QV	1	10	50	400 x 350 x 320	0.045	5.20	1,800
93-773Q	1	10	50	400 x 350 x 320	0.045	4.57	1,800
9-773Q	1	10	50	400 x 350 x 320	0.045	4.57	1,800
9-773QV	1	10	50	400 x 350 x 320	0.045	4.57	1,800

Packaging material	
Primary package	Poly Bag - Low density polyethylene / Insert - 80g Paper
Secondary package	Poly Bag - Low density polyethylene / Sticker - Coated Paper
Tertiary package	Sticker - Coated Paper / Carton - Corrugated Board

Tracheobronchial Nebulizer

Nebulizer and Kits

SINGLE - USE



MEDICAL EQUIPMENT

Storage and handling	
	THIS WAY UP
	Keep away from rain
	Fragile; handle with care
	Keep away from sunlight
	Do not re-use
	Do not use if package is damaged and consult electronic instructions for use
	Consult instructions for use or use electronic instructions for use

SPONSORED IN AUSTRALIA BY

Unoquip Pty Ltd
PO Box 1117
Windsor VIC 3181

1300 55 80 79
info@unoquip.au

DISTRIBUTED BY

UNOQUIP GMBH
Schützengraben 20
8200 Schaffhausen
Switzerland
info@unoquip.com
www.unoquip.com



Ningbo Shengyurui Medical
Appliances Co., Ltd.
No. 138 Binhaisi Rd.
Hangzhou Bay New Zone,
315336, Ningbo City
Zhejiang Province, P.R. China



Shanghai Intl. Holding
Corp. GmbH
Eiffestrasse 80
20537 Hamburg
Germany



MedNet SWISS GmbH,
Bäderstrasse 18,
5400 Baden,
Switzerland